



EMERGENCY TREATMENT COT (ETC) SETUP INSTRUCTIONS

Thank you for purchasing the Kamp-Rite Emergency Treatment Cot
We trust this product will fulfil your absolute requirements in any situation where it is intended for use. Please keep this instruction sheet in a safe place for future setup reference.




1. Remove the ETC from its Carry bag



2. Unfold and position tension straps to inner legs. tighten straps



3. By manipulating the adjustable screws and the tension straps, you can set the ETC up into the desired positions for patient comfort. See Below for setup options.



visit www.kamprite.com for more information on Kamp-Rite Gear